

**EXPERIENCE  
VIVACE.  
EXPERIENCE  
THE RESULTS.™**

**VIVACE**  
Fractional Micro Needle RF

You've just enjoyed  
The Vivace Experience™ ...  
here are tips to optimize your  
recovery and maximize your  
results. No other company is  
more committed to your  
comfort during and after the  
treatment. The Vivace creates  
microchannels that are open  
for four to six hours, creating  
the opportunity to add  
beneficial products during  
and post treatment.\*\*

tip

*Do not reapply makeup right away. If you can, refrain from applying makeup until the following day. Mineral-based makeup is the best post treatment.*

tip

*Do not use any additional products during the first four to six hours following your treatment, unless advised by your clinician.*

STEP

1

**SOOTHE RECOVERY SPRAY:** Immediately following treatment, your clinician will spritz your face with Soothe Recovery Spray. You will take this home with you and should continue to use for the rest of the day as needed.

STEP

2

**ENHANCE WEEKLY MASK:** Your clinician will place the Enhance Weekly Mask on your face immediately following the treatment as part of The Vivace Experience or advise you to use the Enhance Weekly Mask for 20 to 30 minutes within the first four to six hours of the procedure and daily thereafter as directed. To use the mask, remove mask from pouch, unfold mask and apply to your face starting from your nose or center of your face and smooth outwards until your skin is completely covered, shifting the openings for your mouth and eyes to your mouth and eyes. Sit back and Relax.

**SOOTHE RECOVERY SYSTEM:** During the first four to six hours following treatment and for the remainder of the day and evening, use our Soothe Recovery System to support accelerated healing.

STEP  
3

At home continue with Soothe Recovery Spray as outlined in Step 1 above. Spritz your face throughout the day with Soothe Recovery Spray to ease tightness of the face and reduce any remaining and minimal discomfort.

STEP  
4

Depending on the level of your treatment, you may be advised to use Soothe Recovery Balm, which would be applied as needed during the first four to six hours following treatment after you use the Soothe Recovery Spray, and after you have used the Enhance Weekly Mask.

STEP  
5

Before bed, cleanse with Soothe Recovery Cleanser, which contains Azulene, known for its anti-inflammatory benefits. Use very light pressure when cleansing the skin. The cleanser will “melt” off any residue left by the Soothe Recovery Balm.

STEP  
6

After you cleanse with Soothe Recovery Cleanser, you may spritz your face again with Soothe Recovery Spray and use Soothe Recovery Balm as needed.

STEP  
7

You may resume your normal skin routine three days post treatment. We recommend you wait to apply your makeup until the next day.

tip

*Do not exercise or increase your heart rate for 24 hours post treatment, and if instructed by your clinician abstain for up to three days.*

tip

*Avoid all sun exposure.*

tip

*Refrain from deep exfoliation and/or retinol products for seven to 10 days post treatment.*

STEP  
8

We would be honored if you would share your experience with us by providing a testimonial to your physician.

STEP  
9

We recommend you use the Enhance Weekly Mask as part of your regular skincare regime on a weekly basis.

\*\* Your physician or clinician may elect to provide variable instructions for your individual recovery. Please defer to those instructions following the treatment.



THE PERPETUAL PURSUIT OF PERFECTION

AESTHETICS BIOMEDICAL

The innovators at Aesthetics Biomedical are dedicated to providing the best patient experience utilizing the best-in-class device that does it right the first time and delivers results you can trust and share with others.

**[VivaceExperience.com](http://VivaceExperience.com)**